


# Stereotactic Drill System

Product information



Adjustable  
depth stop

Adjustable drill depth

Depth precision: 1 mm

Vibration free,  
low noise

Operated via  
foot switch

Optimised for Deep Brain  
Stimulation (DBS)

Minimally invasive:  
ø 8 mm drill

Minimum tissue warming

Perfect compatibility with  
Precisis Targeting Systems

## Advantages

- Optimal compromise between **minimal** invasion and **sufficient** work space.
- Drill with manual **drill feed** & adjustable **depth stop**.
- Optimise your **functional** stereotactic treatment.
- Perform parallel & coaxial drilling.
- Combine Drill System with our **Target Systems**.

# Stereotactic Drill System

## Product information



### SP-04000 Stereotactic Drill Set

#### Set includes

Product No	Description	Product No	Description
SP-04003-04	Motor mounting tool	SP-04010-A	Foot switch
SP-04007	Drill $\varnothing$ 7 mm	SP-04010-E	Potential equ. cable
SP-04008	Drill $\varnothing$ 8 mm	SP-04011-A	Connection cable
SP-04007-04	Marker Pin	SP-04012	Battery Charger Set
SP-04009-A	Motor unit	SP-04013	Battery, 14.4V, 1.4 Ah AC


#### Compatible Fixation & Targeting Systems



##### CAAT-CSS:

##### Complete Stereotactic System Set

- Center of Arc design
- Cost efficient



##### SP-03000:

##### SP Fixation System Set, CT


- Superior precision Fixation
- Compatible with TA & CAAT






##### TA-07000:

##### TA Target System Set

- Superior precision Targeting
- Analog RM coordinates



#### The System in use



## Summary

The **all new** mechanically guided drill system:  
**Minimally invasive** and optimised for highest accuracy in  
trepanation.  
The perfect tool even for **deep brain stimulation** procedures.

**Maximise**  
**safety and**  
**precision.**